Chemistry Lab Report

Exp No: 01

Exp Name: Standardization of HCl using 100mL of 0.1(M) Na_2CO_3 .

Equipments: 1. Burette(500mL) 2. Conical flask(250mL) 3. Funnel 4. Wash bottle

5. Pipette 6. Burette stand.

Theory: Titration is a common laboratory method of quantitative chemical analysis to determine the concentration of an identified analyte (a substance to be analyzed). A reagent, termed the titrant or titrator, is prepared as a standard solution of known concentration and volume. The titrant reacts with a solution of analyte to determine the analyte's concentration. The volume of titrant that reacted with the analyte is termed the titration volume.

Chemicals: 1. HCl 2. Na_2CO_3 3. Distilled water 4. Indicator(Methyl orange).

Chemical Reaction:

$$\begin{aligned} Na_2CO_3 + 2HCL &\rightarrow 2NaCl + H_2O + CO_2 \\ \text{Now,} \\ e_1S_1V_1 &= e_2S_2V_2 \implies S_2 = \frac{e_1S_1V_1}{e_2V_2} \dots \dots (1) \\ \text{Where, } e_1 &= 2, \quad e_2 &= 1 \\ S_1 &= \text{concentration of } Na_2CO_3 \\ S_2 &= \text{concentration of } HCl \\ V_1 &= \text{volume of } Na_2CO_3 \\ V_2 &= \text{volume of } HCl \end{aligned}$$

Description:

- **01.** $0.1(M)Na_2CO3$ solution preparation: At first we weight 1.061g, then we take it into a 250mL flask and mixed it with 250ml distilled water. After mixing it for some time, it's ready for experimental work.
- **02.** Standardization of HCl solution: At first fil the burette with HCl solution and record the initial burette reading. Then take $10mLNa_2CO3$ in conical flask and mix it up with 1/2 drops of Methyl Orange.
- **03.** Determination: Mix HCl with Na_2CO_3 drop by drop carefully until Na_2CO_3 solution changes it's color. When Na_2CO_3 changes it's color that means it's the end point of our titration. Now mark the reading for HCl from burette for our further calculation.

Data Table:

SN 01

02

10mL

		Burette HCl volume		
V	$V_{1Na_2CO_3}$	$Initial_{(mL)}$	$Final_{(mL)}$	$Diff_{initial-final}$
	10mL	0mL	21.6mL	21.6mL

Calculation:

From the table, Mean value of $HCl=\frac{21.6+21.3}{2}=21.45mL$ From equation (1) we get, $S_2=\frac{e_1S_1V_1}{e_2V_2}=\frac{2\times0.100094\times10}{21.45}=0.09332$ M

49.9mL

21.3mL

Result: The concentration of HCl is **0.09332** M.

21.6mL

Discussion: The molarity of HCl may not be totally correct for same certain chemical fault.

Precaution:

- **01.** Usually an air bubble is present in the nozzele of the burette. It must be removed before taking the inital readint.
- **02.** There should not be any kind of leakage from the burette during titration.
- 03. Always add acid to water.
- **04.** Dont't let base level in burrete to reach zero.