

Chemistry Lab Report

Exp No: 04

Exp Name: Determination of Fe in supplied $FeSO_4$ by using standardized $0.02(M)KMnO_4$ solution.

Equipments: 1. Burette($500mL$) 2. Conical flask($250mL$) 3. Funnel 4. Wash bottle
5. Pipette 6. Burette stand.

Theory: Titration is a common laboratory method of quantitative chemical analysis to determine the concentration of an identified analyte (a substance to be analyzed). A reagent, termed the titrant or titrator, is prepared as a standard solution of known concentration and volume. The titrant reacts with a solution of analyte to determine the analyte's concentration. The volume of titrant that reacted with the analyte is termed the titration volume.

Chemicals: 1. $KMnO_4$ 2. $FeSO_4$ 3. Distilled water 4. Indicator(Methyl orange).

Chemical Reaction:



Description:

- 01.** $0.1(M)KMnO_4$ solution preparation: We obtained $0.02(M)KMnO_4$ provided by our lab assistant.
- 02.** Standardization of $FeSO_4$ solution: At first fill the burette with $KMnO_4$ solution and record the initial burette reading. Then take $10mLKMnO_4$ in conical flask and mix it up with $1/2$ drops of methyl orange.
- 03.** Determination: Mix $KMnO_4$ with $FeSO_4$ drop by drop carefully until $KMnO_4$ solution changes its color. When $FeSO_4$ changes its color that means it's the end point of our titration. Now mark the reading for $KMnO_4$ from burette for our further calculation.

Data Table:

Burette $KMnO_4$ volume				
SN	$FeSO_{4mL}$	$Initial_{(mL)}$	$Final_{(mL)}$	$Diff_{initial-final}$
01	$10mL$	$0mL$	$3.3mL$	$3.3mL$
02	$10mL$	$3.3mL$	$6.6mL$	$3.3mL$
02	$10mL$	$6.6mL$	$10.1mL$	$3.5mL$

Calculation:

From the table, Mean value of $KMnO_4 = \frac{3.3+3.3+3.5}{3} = 3.3666mL$

From equation (1) we get,

$$1000mL \ 1(M) \ KMnO_4 \equiv 5mol \ Fe^{2+}$$

$$1mL \ 1(M) \ KMnO_4 \equiv \frac{5}{1000}mol \ Fe^{2+}$$

$$3.36mL \ 0.02(M) \ KMnO_4 \equiv \frac{5 \times 3.36 \times 0.02 \times 55.85}{1000}g \ Fe^{2+}$$

Therefore,

$$10mL \ Fe^{2+} = 0.0187656g \ Fe^{2+}$$

$$100mL \ Fe^{2+} = \frac{0.0187656 \times 100}{10} = 0.187656g \ Fe^{2+}$$

Result: The amount of Fe^{2+} is **0.0.187656g**.

Discussion: The amount of Fe may not be totally correct for same certain chemical fault.

Precaution:

- 01.** Usually an air bubble is present in the nozzle of the burette. It must be removed before taking the initial reading.
- 02.** There should not be any kind of leakage from the burette during titration.
- 03.** Always add acid to water.
- 04.** Don't let base level in burette to reach zero.